



PRO SIMBODIES - Thoracotomy Torso

PRODUCT NUMBER: 1383



DESCRIPTION:

PRO SIMBODIES **Thoracotomy Trainer** is a **procedural, repairable and operable biofidelic manikin** that provides new immersive training options for clinical educators and practitioners of emergency medical care.

Adaptable, portable and realistic, this manikin can deliver surgical training opportunities in a range of austere **emergency training environments**. It can be cut open, operated on and sutured back together using standard surgical equipment.

SKILLS:

- To educate and realistically practice clamshell thoracotomy, pericardial tamponade and suturing
- Needle Decompression

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- Chest Drain
- Available in both a full body and torso task trainer
- Closed eye lids capable of manual opening
- Ability to perform basic airway training

CHARACTERISTICS:

GENERAL:

- Anatomically correct airways, and oesophagus, which can be intubated and lungs that can be ventilated.
- Permanent wound to the heart wall with manual bleeding capabilities
- The beating heart (manually pumped) and pericardium can simulate a realistic cardiac tamponade and can be accessed via a fully repairable clamshell thoracotomy.
- Feature a realistic thoracic skeleton a replaceable sternum and detailed skin, fat and subcutaneous tissues, all complete with varying palpability and realistic colour.
- It can be cut open, operated on and sutured back together using standard surgical equipment.
- Replaceable and repairable features including replaceable skins and pericardium

TECHNICAL SPECIFICATIONS:

- Splashproof and washable skin
- Robust & re-usable
- Skin safe, Prosthetic Grade Platinum Silicone
- Age: Approx. 30-40 years old
- Male only

INCLUDE:

- 1 Pro Simbodies Torso

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>