



FREDDIE FISTULA

PRODUCT NUMBER: 1453



DESCRIPTION

Freddie Fistula Skills Trainer was developed to provide clinicians with a realistic aid to facilitate instruction, practice, and management of open abdominal wounds with active fistulas.

The Freddie Fistula Skills Trainer incorporates a realistic feeling soft abdomen with an open abdominal wound showing exposed intestines with "gutters" on all sides, a stoma and also a drainage tube.

SKILLS

Learn multiple wound care modalities to control fistula effluent and allow for optimal wound healing.
Build competency and confidence in this challenging wound care scenario
Fistula dressing scenarios

CHARACTERISTICS

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 671 222 375 / (+351) 917 260 661
www.morethansimulators.com



Three separate fistulas are present in the abdominal wound area:

Stomatized Fistula ? Easy to observe effluent flow

Non-Stomatized Fistula ? Difficult to observe when effluent is not flowing

Hidden Fistula ? Location is not observed

Any combination of the three fistulas can permit effluent flow via color-coded pinch clamps, allowing for a challenging scenario when applying a dressing.

INCLUDE

1 Freddie fistula skills trainer

1 Base with fluid reservoir

1 Carrying bag (option)